

Athinoula A.

Martinos
Center
For Biomedical Imaging



MASSACHUSETTS
GENERAL HOSPITAL



HARVARD
MEDICAL SCHOOL



Massachusetts
Institute of
Technology

Multiparametric mapping with Pulseq

Shohei Fujita

BRAIN (Bilgic Reconstruction Acquisition for Imaging Neuroscience) Lab

Objectives

1. Demonstrate the actual workflow of implementing and executing custom pulse sequences in Pulseq
2. Introduce an application of Pulseq for cross-platform, multiparametric mapping

Background


- Multi-site data cannot be naively pooled together for analysis due to large inter-scanner differences

Address the unmet need for a cross-platform, multiparametric technique to facilitate data harmonization across different sites

3D-QALAS: Rapid whole-brain T1 and T2 mapping

3D-QALAS^{1,2}:

- 3D-quantification using an interleaved Look–Locker acquisition sequence with a T2 preparation pulse
- Multi-acquisition 3D gradient echo
- Signal modeling-based parameter fitting on the five source images to derive quantitative maps




[1] Kvernby et al. *J Cardiovasc Magn Reson* 2014; [2] Fujita et al. *Magn Reson Imaging* 2019

3D-QALAS: Rapid whole-brain T1 and T2 mapping


3D-QALAS^{1,2}:

- 3D-quantification using an interleaved Look–Locker acquisition sequence with a T2 preparation pulse
- Multi-acquisition 3D gradient echo
- Signal modeling-based parameter fitting on the five source images to derive quantitative maps



[1] Kvernby et al. *J Cardiovasc Magn Reson* 2014; [2] Fujita et al. *Magn Reson Imaging* 2019


Study overview




Site	Manufacturer	Scanner model	Software version	Evaluation
Site A	Siemens Healthcare, Erlangen, Germany	Prisma	VE11C	-
Site A	Siemens Healthcare, Erlangen, Germany	Prisma	XA30	Cross-software
Site A	Siemens Healthcare, Erlangen, Germany	Skyra	VE11C	Cross-scanner
Site B	Siemens Healthcare, Erlangen, Germany	Prisma	VE11C	Cross-site
Site C	GE HealthCare, WI, USA	Premier XT	MR30	Cross-vendor

Sequence implementation on Pulseq

writeSeq.m



.seq




```
# Pulseq sequence file
# Created by writeSeq.m toolbox
[VERSION]
major 4
minor 4
revision 0

[DEFINITIONS]
AddRasterTime 1e-05
BlockDurationInter 1e-05
FOV 0.1502 0.1502 0.1502
GradientInterpolate 0.45
RadiofrequencyRasterTime 1e-05
TotalDuration 261

# Format of blocks:
# MAM DUR RF GX GY GZ ADC EXT
[BLOCKS]
1 1000 1 0 0 0 0 0 0
2 2500 2 0 0 0 0 0 0
3 2500 2 0 0 0 0 0 0
4 2500 2 0 0 0 0 0 0
5 1500 2 0 0 0 0 0 0
6 42 3 0 0 0 0 0 0
7 970 0 0 0 1 0 0
8 32 4 0 0 0 0 0 0
9 10 5 0 0 0 0 0 0
10 58 0 2 3 0 0 0 0
11 306 0 5 0 0 0 1 0
12 180 0 5 0 0 0 0 0
13 10 0 0 0 0 0 0 0
14 22 5 0 0 0 0 0 0
15 7 0 0 0 0 0 0 0
16 10 0 0 0 0 0 0 0
17 306 0 5 0 0 2 0 0
18 180 0 6 18 0 0 0 0
19 7 0 0 0 0 0 0 0
20 2 0 0 0 0 0 0 0
21 7 0 0 0 0 0 0 0
22 58 0 2 11 0 0 0 0
23 306 0 5 0 0 3 0 0
24 10 0 0 0 12 0 0 0
25 7 0 0 0 0 0 0 0
26 22 7 0 0 0 0 0 0
27 7 0 0 0 0 0 0 0
28 58 0 2 13 0 0 0 0
29 306 0 5 0 0 4 0 0
```

Scanner



Sequence implementation on Pulseq: writeSeq.m

```
%-----  
% Define high-level parameters  
%-----  
  
% Set system parameters  
sys = mr.opts('MaxGrad', 32, 'GradUnit', 'mT/m', ...  
    'MaxSlew', 160, 'SlewUnit', 'T/m/s', ...  
    'rfDeadTime', 100e-6, ...  
    'rfRingdownTime', 60e-6, ...  
    'adcDeadTime', 40e-6, ...  
    'adcRasterTime', 2e-6, ...  
    'gradRasterTime', 20e-6, ...  
    'blockDurationRaster', 20e-6, ...  
    'B0', 3.0); |  
  
gradRasterTime (Siemens) = 10us  
gradRasterTime (GE) = 4us  
  
% Create sequence object and define FOV, resolution, and other high-level  
seq = mr.Sequence(sys); % Create a new sequence object  
fov=[192 168 168]*1e-3; % Define FOV and resolution  
N = [192 168 56];  
  
Nx = N(1);  
Ny = N(2);  
Nz = N(3);  
  
bandwidth = 340; % Hz/pixel  
flip_angle = 4; % degrees  
rfSpoilingInc=117; % RF spoiling increment  
dwell = 16e-6;  
Tread = dwell*Nx;
```

Sequence implementation on Pulseq: writeSeq.m

```
%-----  
% Prepare sequence blocks  
%-----  
  
deltak = 1 ./ fov;  
gx = mr.makeTrapezoid('x',sys,'Amplitude',Nx*deltak(1)/Tread,'FlatTime',ceil(Tread/sys.gradRasterTime)*sys.g  
gxPre = mr.makeTrapezoid('x',sys,'Area',-gx.area/2); % Gx prewinder  
gxSpoil = mr.makeTrapezoid('x',sys,'Area',gx.area); % Gx spoiler  
  
% Make trapezoids for inner loop to save computation  
  
gyPre = mr.makeTrapezoid('y','Area',deltak(2)*(Ny/2),'Duration',mr.calcDuration(gxPre), 'system',sys);  
gzPre = mr.makeTrapezoid('z','Area',deltak(3)*(Nz/2),'Duration',mr.calcDuration(gxPre), 'system',sys);  
  
gyReph = mr.makeTrapezoid('y','Area',deltak(2)*(Ny/2),'Duration',mr.calcDuration(gxSpoil), 'system',sys);  
gzReph = mr.makeTrapezoid('z','Area',deltak(3)*(Nz/2),'Duration',mr.calcDuration(gxSpoil), 'system',sys);  
  
stepsY=((traj(:,2)-1)-Ny/2)/Ny*2;  
stepsZ=((traj(:,1)-1)-Nz/2)/Nz*2;  
  
% Create non-selective pulse  
rf = mr.makeBlockPulse(flip_angle*pi/180, sys, 'Duration', 0.1e-3);
```


Sequence implementation on Pulseq: writeSeq.m

```
% Prep pulses
% T2 prep and IR prep pulse are imported from external txt files
% txt file is in mag and phase, while mr.makeArbitraryRf assumes real and
t2prep = readmatrix('csv_imports/T2prep.txt');
Re = t2prep(:,1) .* cos(t2prep(:,2));
Im = t2prep(:,1) .* sin(t2prep(:,2));
t2prep_pulse = mr.makeArbitraryRf((Re+Im*1i).', 380.4*pi/180, 'system',sy


IRprep = readmatrix('csv_imports/IRprep.txt');
Re = IRprep(:,1) .* cos(IRprep(:,2));
Im = IRprep(:,1) .* sin(IRprep(:,2));
IRprep_pulse = mr.makeArbitraryRf((Re+Im*1i).', 850*pi/180, 'system',sys,
```

mr.makeArbitraryRF

T2prep



IR



Sequence implementation on Pulseq: writeSeq.m

```
%-----  
% Build sequence  
%-----  
  
for iZ = 1:nTR  
    rf_phase=0;  
    rf_inc=0;  
  
    for iY = 1:nETL  
  
        % Calculate index for ky-kz look up table  
        index = iY + (iZ-1)*nETL;  
  
        % RF spoiling  
        rf.phaseOffset=rf_phase/180*pi;  
        adc.phaseOffset=rf_phase/180*pi;  
        rf_inc=mod(rf_inc+rfSpoilingInc, 360.0);  
        rf_phase=mod(rf_phase+rf_inc, 360.0);      %increment RF phase  
  
        % Excitation  
        if iY == 1  
            seq.addBlock(rf,mr.makeLabel('SET', 'LIN', segmentID)); % segment  
        else  
            seq.addBlock(rf);  
        end  
  
        % Encoding  
        seq.addBlock(gxPre, ...  
            mr.scaleGrad(gyPre,-stepsY(index)), ...  
            mr.scaleGrad(gzPre,-stepsZ(index)));    % Gz, Gy blips, Gx pre-win
```

Sequence implementation on Pulseq: timing check, PNS check

```
%-----
% Check timing and write sequence
%-----

% check whether the timing of the sequence is correct
[ok, error_report]=seq.checkTiming;

if (ok)
    fprintf('Timing check passed successfully\n');
else
    fprintf('Timing check failed! Error listing follows:\n');
    fprintf([error_report{:}]);
    fprintf('\n');
end

% PNS check
[pns_ok, pns_n, pns_c, tpns]=seq.calcPNS('/autofs/cluster/berkin/fujita/mapping/seq');



if (pns_ok)
    fprintf('PNS check passed successfully\n');
else
    fprintf('PNS check failed! The sequence will probably be stopped by timing\n');
end

% Set definitions for reconstruction
seq.setDefinition('FOV', fov);
seq.setDefinition('Matrix', N);
seq.setDefinition('nETL', nETL);
seq.setDefinition('nTR', nTR);
seq.setDefinition('traj_z', traj(:,1));
seq.setDefinition('traj_y', traj(:,2));

% plot
seq.plot('TimeRange',[0 1.5], 'timeDisp', 'ms');

% Write to pulseq file
filename = strrep(mfilename, 'write', '');
seq.write([filename,'.seq']);
```

Sequence implementation on Pulseq



FOV: 192x168x168

Matrix: 192x168x56

TR: 4500ms

TE: 2.3ms


Echo spacing: 5.8ms

FA: 4°


ETL: 127

TA: 4:21


Center-out elliptical ordering
Cartesian sampling




Before going to the scanner... simulating the .seq file in POET



Running your sequence on the scanner



Running your sequence on the scanner



Exporting and reconstructing your pulseq data

Siemens

Twix



meas_MID*.dat

```
data_file_path =  
'/autofs/cluster/berkin/fujita/testscan/archive/meas  
_MID00465_FID38558_pulseq_qalas.dat';
```

```
twix_obj = mapVBVD2(data_file_path);  
data_unsorted = twix_obj{end}.image.unsorted();
```

```
% size(data_unsorted): adc_len, ncoil, readouts
```



Data sorted in the order it was acquired
Sort using your seq file

GE


```
pfile = 'P06144.7';
```

```
data_unsorted =  
toppe.utils.loadpfile(pfile, [], [], [], 'acqOrder', true);
```

```
% size(data_unsorted): adc_len, ncoil, readouts
```




Cross-platform multiparametric mapping








Identical implementation and reconstruction pipeline improves cross-vendor reproducibility


(A) Cross-vendor reproducibility




(B) Comparison of Pulseq- and vendor-native implementation




Implementation for in vivo acquisition




FOV: 256x240x208 (ADNI-compliant)
Matrix: 256x240x208
TR: 4500ms
TE: 2.3ms
Echo spacing: 5.8ms




Vendor-agnostic implementation improves reproducibility



(B) Brain structure T1 value reproducibility




(C) Brain structure T2 value reproducibility



Summary

- Workflow for using custom sequences implemented in Pulseq was demonstrated
- An application of multiparametric mapping using Pulseq was demonstrated
- Pulseq implementation exhibited:
 - Higher reproducibility than vendor-native implementations
 - ADNI-compliant field-of-view sizes with 1mm isotropic resolution within 5 minutes, while maintaining a cross-platform coefficient of variation below 4%.



Athinoula A.

Martinos Center

For Biomedical Imaging



MASSACHUSETTS
GENERAL HOSPITAL



HARVARD
MEDICAL SCHOOL



Massachusetts
Institute of
Technology

Thank you!

Shohei Fujita

BRAIN (Bilgic Reconstruction Acquisition for Imaging Neuroscience) Lab